

Attune® NxT 3激光11色荧光选择表

Excitation Laser	Emission Filter (nm)	Recommended Dyes	Fluorescent Proteins
Violet-405 nm	VL1: 440/50	Alexa Fluor 405, eFluor 450, SB436, Pacific Blue, V450, BV421	eCFP, eBFP
	VL2: 512/25	eFluor 506, Pacific Green, V500, BV510	vGFP
	VL3: 603/48	SB600, Pacific Orange, Qdot 605, BV605, eFluor 605NC, eFluor 625NC	
	VL4: 710/50	SB702, Qdot 705, BV711	
Blue-488 nm	BL1: 530/30	Alexa Fluor 488, FITC, CFSE	eGFP, eYFP
	BL2: 574/26	PE, PI	
	BL3: 695/40	PerCP-eFluor 710, PerCP-Cy5.5, PE-Cy5.5, PE-Cy5, 7-AAD, Qdot 705, PI, PerCP, PE-Alexa Fluor 700	
	BL4: 780/60	PE-Cy7, Qdot 800	
Red-637 nm	RL1: 670/14	APC, Alexa Fluor 647, eFluor 660, Qdot 655	
	RL2: 720/30	Alexa Fluor 700, Alexa Fluor 680, APC-Alexa Fluor 700, Qdot 705, APC-R700	
	RL3: 780/60	APC-Alexa Fluor 750, APC-eFluor 780, APC-Cy7, Qdot 800, APC-H7	

